ic1011

THERAPEUTIC MASSAGE CHAIR

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USER MANUAL

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**PLEASE DO NOT RETURN
THIS PRODUCT TO THE STORE**

**FOR QUESTIONS OR ASSISTANCE ON ASSEMBLY, PARTS OR OPERATION,
CONTACT OUR CUSTOMER SERVICE!**

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1 866 967.7333 ext. 227

Thank you for purchasing this product. Please read this instruction manual carefully before operating this massage chair.

Please pay attention to the important safety warnings and keep this manual for future use.

Note: EBI reserves the right to revise the design and description of this product without notice.

1. CONTENTS

1. Safety warnings	2
2. Environment.....	2
3. Restrictions.....	2
4. Maintenance.....	3
5. Troubleshooting.....	3
6. Features and benefits	3
7. Installation.....	4
8. Components	7
9. Operation	9
10. Floor protection.....	13
11. Product specifications.....	13
12. Warranty	14
FRANÇAIS	15

1. SAFETY WARNINGS

- This product is not recommended for children.
- Disconnect the chair before cleaning.
- Please pull out the plug after use or before cleaning to avoid injury or damage to the product.
- Please operate this product according to the instructions in this manual.
- Use only accessories recommended by iComfort®.
- This product is designed for indoor use only.
- Recommended use time is 30 minutes or less.
- Please avoid use if the chair is damaged and call your iComfort® technician or dealer.
- Make sure the chair radiator in the lower back behind the chair is well ventilated and not obstructed.
- Please make sure small objects are not dropped into the chair mechanism.
- Use of this chair under the influence of drugs or alcohol is not recommended.
- Do not use wet hands to pull the AC cord.
- Do not spill liquids onto this product to avoid electric shock or cause damage to this product.
- Maintain your chair and avoid damaging the wires.
- Do not use a wet cloth to clean the electric parts.
- In the event of a power failure disconnect the chair.
- If the chair works abnormally, stop using it, disconnect the chair and call your iComfort® service technician or customer service representative.

2. ENVIRONMENT

This chair should be used only in a well ventilated, normally dry and dust free environment.

3. RESTRICTIONS

People with osteoporosis, heart disease and heart pacemakers, pregnancies, injuries or any other health related issues should consult their doctor before using this massage chair.

Stop immediately the use of this massage chair if you feel unwell and consult your doctor.

4. MAINTENANCE

- This product should only be maintained by iComfort® certified technicians. Otherwise, EBI reserves the right to void the warranty.
- Disconnect the power cord after use.
- Make sure the wall outlet is in good working order before use.
- If this chair will not be used for an extended period of time, store in a dry and dustless environment.
- Avoid chair exposure to high temperatures and direct sunlight.
- Clean this product with a dry cloth. Do not use thinner, benzene or alcohol.
- When moving the chair on wood floor, rough ground or in a narrow space two people should lift the chair together.
- We recommend to use the massage chair intermittently and not for long time periods.

5. TROUBLESHOOTING

- It is normal for the motor to make noise while in use.
- If the controller doesn't work normally, please check that the plug and socket are properly connected and the switch is turned on.
- When the programmed working time is over, the chair will power-off automatically.
- After several use, the chair might power-off because of overheat. It can be used again after half an hour of cooling time.

6. FEATURES AND BENEFITS

Enjoy a full en invigorating massage with the iComfort® ic1011 massage chair. You will feel as if you had your very own massage therapist. The Body Scanning Technology actually maps the back of each individual user in order to provide a customized massage. The handy retractable footrest gives you maximum configuration options. The iComfort® ic1011 massage chair reproduces over 5 massage actions all with the sensation of a professional massage in the comfort of your very own home! The iComfort® ic1011 massage chair will provide you with many years of soothing massage and comfort.

- One set of vertically movable quad-rollers assembled with high quality soundproofed mechanical massage hands.
- Five massage functions: Shiatsu kneading, flapping, knocking, pressing and simultaneous kneading flapping.
- Automatic spinal cord scans to determine the best position for the massage.
- 11 air bags for the lower body are conducted by with four different modes and two intensity levels.
- LCD remote control shows timers, massage functions, action modes, massage speed, intensity, width of the massage hands and massage positions.
- The massage hands will reposition automatically after the massage stops.
- The leg massage section can be power-lifted and reclined independently or combined with the backrest.

7. INSTALLATION

The ic1011 massage chair comes partly assembled.


- Assemble this massage chair according to the installation instructions provided in this guide.
- Check the parts list carefully to make sure nothing is missing.
- Sitting in the massage chair before it is completely assembled could result in permanent damage to the chair and personal injury.
- **This massage chair is equipped with a cord having a grounding conductor. The cord must be plugged into an appropriate AC wall outlet.**
- If the massage chair fuse is burned out, please turn the chair off and unplug the power. You can unscrew the fuse located in the switch box at the bottom left of the chair in the back.
 - A) Unplug the AC connector of the chair from the wall outlet.
 - B) Gently unscrew the fuse cover and the fuse will come out at the same time

Use a same size and type of fuse to replace the old one and then cover the lid. You can call iComfort® at 1-866-967-7333 ext 227 if you need help.

CAUTION: Never substitute a missing part with one of your own. If a part is missing, contact customer service at: **1-866-967-7333 ext 227**




Unpacking

1. Remove the massage chair, cushions and accessories from the main box.
2. Remove all plastic wraps.
3. Slowly lift the backrest up to approximately 120 degrees until you hear a light 'click' sound. This signifies the backrest is secured to the seat area.

Remote control and holder


Install remote control and holder. Attach the remote control holder underneath the right armrest with the screws provided.

Seat air hoses

There are three air connectors under the seat cushion. Connect them to the air connectors on the massage chair seat. Please make sure to connect them using the same matching codes. Once done, put the seat back in.


IMPORTANT: Please make sure not to bend the airbag hoses. This could prevent air from getting to the airbags and cause a malfunction.

Back and head cushions

Install the back and head cushions using the cushion's Velcro tape.

The back and its cushion are connected by a zipper

The pillow pad and the back cushion are connected by a velcro strip


Final Installation

1. Insert the AC power cord into the wall outlet
2. Turn on the main power switch at the bottom and behind the chair on the control box.
3. Use the remote control to choose your preferred massage function and enjoy the comfortable massage.

8. COMPONENTS

Vent


Power

Fuse box

Wire

Plug

Caster


Massage Rollers

Massage
mechanism
movement

Seat air bags

Footrest
airbags

9. OPERATION

First, adjust the massage chair angle. Then, press on the power button to start the massage. The ic1011 will scan the length of your back and width of your shoulders.

We recommend adjusting the chair angle before powering up the chair. This way, the body measurements will be done correctly. Once done (about 20 seconds), the automatic massage mode (mode A) will start. You can press on the 'Mode button repeatedly to choose another automatic massage mode. There are four modes (A,B,C and D). These massage modes will go through most of the massage techniques available with this chair. After the massage is done, the chair will return to its original position. You can interrupt the massage at any point by shutting off the unit. The massage chair will then return to its original position.

If you prefer, you can also adjust a preferred shoulder positions. After pressing on the power on button, press on the 'Adjusting shoulder first' button. There are five choices (from low to high). You can also choose a massage area. The automatic mode will then concentrate on this area. Press one of the three buttons to choose which area of your back you want the massage done (Upper body, Neck shoulder or Waist back). The massage chair will go through its multiple massage modes for about 15 minutes.

If you prefer a specific massage mode, you need to use the manual massage mode. If so, you can choose the type of massage you want (knead, flap, etc.). The massage chair will then use the massage mode you have chosen.

There are other manual massage modes like the 'Arthrosis traction' (knee stretching), footrest airbag massage and seat massage. All these can be activated while in automatic or manual mode.

You can also adjust massage speed, massage mechanism width and airbag intensity at any time.

Button Functions

Press the arrow up button to raise the footrest only. Press the arrow down to recline the foot rest section only..


Press the arrow up button to raise the footrest only. Press the arrow down to recline the foot rest section only. Press the falling or rising button to adjust the massage chair angle. This feature raises the footrest and reclines the back section at the same time.

Press on the power button to start the massage. Press again and the massage chair will come back to its original position.

Press on the 'adjusting shoulder first' to adjust the chair to your preference. Press on the arrows (up or down).


Press one of the three buttons to choose which area of your back you want the massage done (Upper body, Neck shoulder or Waist back). For Upper body (whole back), the massage device works with various techniques at different body zones or your whole back. For Neck shoulder, the massage device automatically covers you neck or shoulder. For Back waist, the massage device automatically covers your back and waist.


Manual mode

Press this button to activate the Kneading massage feature.

Press this button to activate the Flapping massage feature.


Press this button to activate the Kneading and Flapping massage features at the same time.

Press this button to activate the Pressing massage feature.

Press this button to activate the Knocking massage feature.

Upper body massage adjustments

There are five speed massage options available.

There are three width adjustment options for flapping, kneading or knocking. This feature adjusts the width of the massage mechanism.

Press to select overall, partial or fixed back stretching massage. Press repeatedly to choose.

Press to select back stretching massage. Press the arrow up button to move the massage mechanism up to suit a precise massage points. Press the arrow down to move the massage device downward to suit precise massage points.

Press this button to start air massage for the lower body.

Press this button to choose a full massage mode. There are four modes (A,B,C and D). These massage modes will go through most of the massage features available with this chair.


Press this button to choose different intensities of air massage.

Press this button to activate the arthrosis traction (knee stretching) massage feature.


Moving the chair

Make sure that all wires are away from the ground, turn the back backward to a proper angle with the gravity center of the chair falls onto the casters. Pull the massage chair forward or backward as you wish, and then turn the back to its normal angle.

10. FLOOR PROTECTION

Prolonged use of this massage chair can damage wood floors. In order to prevent permanent damage, place a carpet or mat under the chair.

11. PRODUCT SPECIFICATIONS

Net Weight: 83kg

Dimensions: Reclined: 170 cm length x 80 cm width x 77 cm High

Upright: 111 cm length x 80 cm width x 125 cm High

Power: 180w

Voltage: 110V, 60 Hz

12. WARRANTY

LIMITED ONE (1) YEAR WARRANTY

This product is guaranteed against all manufactured faults and defects for a period of One (1) year from the date of purchase.

For repairs or product defects, please call our customer service department at: **1-866-967-7333**. Please have your proof of purchase on hand. If we are not able to replace your defective part, we will send you a service technician for repairs.

Important notes: Any damages caused by an accident or an abusive usage are **NOT** covered under warranty. We also decline all responsibility concerning indirect or consecutive damages. In certain Canadian provinces, it is not permitted to limit indirect or consecutive damages. Therefore it is possible these limitations are not applicable to you. This warranty gives you specific legal rights. You may also have other legal rights which vary from one province to another.

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You can also fill out the registration card on our website at **www.iComfort.ca**

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Pictures are for reference only, colors and models may vary without notice. Made in China. iComfort™ is a trademark of EBI Inc. Photographies pour référence seulement. Les couleurs et modèles sont sujet à varier sans préavis. Fabriqué en Chine iComfort^{MC} est une marque de commerce de EBI inc